


Surgical sutures for

SKIN CLOSURE

Surgical
sutures



**SERAG
WIESSNER**



Convincing down to the smallest detail

Decades of experience

High-tech and hand-crafted

For nearly 150 years, SERAG-WIESSNER has stood for surgical sutures of the highest quality. As the oldest established German manufacturer, the company combines decades of tradition with state-of-the-art technology. Extremely precise automated production technology is combined with the reliable and diligent work of our experienced employees. In addition, we work closely together with medical professionals in both clinical practice and research. This results in an extensive range of top-quality suture materials for skin closure.

Non-absorbable sutures

SERALON® (blue and undyed)

Compared to other monofilament suture material, this non-absorbable polyamide thread is noticeably soft and smooth. It also possesses the tying effect important for knot security. Thanks to its ease

of use, SERALON® has been successfully used as a skin suture for decades.

SUPRAMID (black and undyed)

SUPRAMID has excellent knot placement and is greatly valued for its lack of capillarity (comparable to monofilaments).

SERAPREN® (blue)

This monofilament, non-absorbable, polypropylene thread can also be recommended for skin sutures due to its smooth, non-porous surface.

SERALENE® (blue)

Due to its smooth, non-porous surface, this non-absorbable thread guarantees optimal sliding properties in the tissues. The substance PVDF is related to PTFE.

Absorbable threads

SERAFIT® (violet and undyed) / SERAPID® (undyed)

These materials are braided, absorbable threads made from polyglycolic acid with a special coating. They are suitable for intracutaneous skin closure.

SERASYNTH® (violet)

This polydioxanone thread belongs to the group of absorbable materials. The monofilament is characterised by its excellent sliding properties and high linear and knot tensile strength.

Intracutaneous suture incutifix®


incutifix® intracutaneous sutures are used where optimal cosmetic effects are required. At both the needle entry and exit points, the non-absorbable, monofilament thread is fixed with a

SERAFAST® (violet and undyed)

This synthetic and absorbable monofilament thread made from polyglycolic acid-caprolactone is particularly suitable for indications with a short wound-healing phase.


Suture material

Intracutaneous suture



incutifix® - a considerable simplification of thread anchoring of intracutaneous sutures through the use of special clips

clip. The thread runs parallel to the surface through the corium.



SERAG-WIESSNER 300 steel needles have very high ductility and low penetrating forces.


Needles

Needle shape
D = $\frac{3}{8}$ circle


Needle cross-section
S = cutting needle

Special characteristics
S = slim point

Example: DSS 18
 $\frac{3}{8}$ circle
cutting needle
slim point
18 mm long (straightened length)



Reference product 1
Reference product 2
SERAG-WIESSNER 300-series steel needles



Needles

Our atraumatic needles are stainless steel needles with high flexural resistance and outstanding penetration capacity, thus permitting safe and easy work.

The designation of our atraumatic needles uses a combination of letters and numbers as recommended by the Technical Committee of the Association of Surgical Suture Manufacturers. The first two letters denote the needle shape and cross-section. A third or even fourth letter indicates special characteristics of the needle.

The number following the letters gives the needle length (true length when straightened out).

Penetration force

DSS needles, our best needle generation for skin closure, have a particularly sharp point. This leads to considerably lower penetration forces and much smaller puncture sites than with conventional needles. The diagrams below left illustrate the corresponding comparisons made by our test laboratory.



reverse cutting needle



reverse cutting needle with special point

Needle designation

DS-10



Material	Size/USP	Size/EP	Length of suture/metre	Units per box	Code number undyed	Code number dyed
SERALON	7/0	0,5	1 x 0,50	24	-	VO051711
SERALON	6/0	0,7	1 x 0,50	24	-	VO071711
SERALON	6/0	0,7	1 x 0,75	24	-	VO073411
SERALON	5/0	1	1 x 0,50	24	-	VO101711
SERALON	5/0	1	1 x 0,75	24	-	VO103411
SERALON	4/0	1,5	1 x 0,50	24	-	VO151711
SERALON	4/0	1,5	1 x 0,75	24	-	VO153411

DS-12



SERAFIT	5/0	1	1 x 0,45	36	50101312	-
SERALENE	5/0	1	1 x 0,50	24	-	LO101712
SERALENE	4/0	1,5	1 x 0,50	24	-	LO151712
SERALON	6/0	0,7	1 x 0,50	24	-	VO071712
SERALON	6/0	0,7	1 x 0,75	24	-	VO073412
SERALON	5/0	1	1 x 0,50	24	-	VO101712
SERALON	5/0	1	1 x 0,75	24	-	VO103412
SERALON	4/0	1,5	1 x 0,50	24	-	VO151712
SERALON	4/0	1,5	1 x 0,75	24	-	VO153412
SERALON	3/0	2	1 x 0,50	24	-	VO201712
SERALON	3/0	2	1 x 0,75	24	-	VO203412
SERAPREN	6/0	0,7	1 x 0,50	24	-	CO071712
SUPRAMID	5/0	1	1 x 0,50	24	-	TO101712
SUPRAMID	4/0	1,5	1 x 0,50	24	-	TO151712

DS-15



SERAFIT	4/0	1,5	1 x 0,70	24	40153213	-
SERAFIT	6/0	0,7	1 x 0,45	36	50071313	-
SERAFIT	5/0	1	1 x 0,45	36	50101313	-
SERAFIT	5/0	1	1 x 0,70	36	-	60103213
SERAFIT	4/0	1,5	1 x 0,45	36	50151313	-
SERALENE	5/0	1	1 x 0,50	24	-	LO101713
SERALENE	5/0	1	1 x 0,75	24	-	LO103413
SERALENE	4/0	1,5	1 x 0,50	24	-	LO151713
SERALENE	3/0	2	1 x 0,50	24	-	LO201713
SERALENE	3/0	2	1 x 0,75	24	-	LO203413
SERALON	6/0	0,7	1 x 0,50	24	-	VO071713
SERALON	6/0	0,7	1 x 0,75	24	-	VO073413
SERALON	5/0	1	1 x 0,35	24	-	VO100613
SERALON	5/0	1	1 x 0,50	24	-	VO101713
SERALON	5/0	1	1 x 0,75	24	-	VO103413
SERALON	4/0	1,5	1 x 0,50	24	-	VO151713
SERALON	4/0	1,5	1 x 0,75	24	-	VO153413
SERALON	3/0	2	1 x 0,50	24	-	VO201713
SERALON	3/0	2	1 x 0,75	24	RO203413	VO203413
SERALON	2/0	3	1 x 0,50	24	-	VO301713
SERALON	2/0	3	1 x 0,75	24	-	VO303413
SERAPID	4/0	1,5	1 x 0,45	36	30151313	-
SERAPREN	5/0	1	1 x 0,50	24	-	CO101713
SERAPREN	4/0	1,5	1 x 0,50	24	-	CO151713

Continued on next page!

Needle designation

DS-18

Material	Size/USP	Size/EP	Length of suture/metre	Units per box	Code number undyed	Code number dyed
SERASYNTH	5/0	1	1 x 0,70	24	-	90103213
SERASYNTH	4/0	1,5	1 x 0,70	24	-	90153213
SUPRAMID	6/0	0,7	1 x 0,50	24	-	TO071713
SUPRAMID	5/0	1	1 x 0,50	24	-	TO101713
SUPRAMID	4/0	1,5	1 x 0,50	24	-	TO151713
SUPRAMID	3/0	2	1 x 0,50	24	-	TO201713
SERAFAST	5/0	1	1 x 0,45	24	40101314	-
SERAFAST	5/0	1	1 x 0,70	24	40103214	JO103214
SERAFAST	4/0	1,5	1 x 0,45	24	40151314	-
SERAFAST	4/0	1,5	1 x 0,70	24	40153214	-
SERAFAST	3/0	2	1 x 0,45	24	40201314	-
SERAFIT	5/0	1	1 x 0,45	36	50101314	-
SERAFIT	4/0	1,5	1 x 0,45	36	50151314	60151314
SERAFIT	4/0	1,5	1 x 0,70	36	50153214	60153214
SERAFIT	3/0	2	1 x 0,45	36	50201314	60201314
SERAFIT	3/0	2	1 x 0,70	36	50203214	60203214
SERAFIT	2/0	3	1 x 0,70	36	50303214	-
SERALENE	5/0	1	1 x 0,50	24	-	LO101714
SERALENE	5/0	1	1 x 0,75	24	-	LO103414
SERALENE	4/0	1,5	1 x 0,50	24	-	LO151714
SERALENE	4/0	1,5	1 x 0,75	24	-	LO153414
SERALENE	3/0	2	1 x 0,50	24	-	LO201714
SERALENE	3/0	2	1 x 0,75	24	-	LO203414
SERALENE	2/0	3	1 x 0,50	24	-	LO301714
SERALON	6/0	0,7	1 x 0,50	24	-	VO071714
SERALON	6/0	0,7	1 x 0,75	24	-	VO073414
SERALON	5/0	1	1 x 0,35	24	-	VO100614
SERALON	5/0	1	1 x 0,50	24	-	VO101714
SERALON	5/0	1	1 x 0,75	24	-	VO103414
SERALON	4/0	1,5	1 x 0,35	24	-	VO150614
SERALON	4/0	1,5	1 x 0,50	24	RO151714	VO151714
SERALON	4/0	1,5	1 x 0,75	24	-	VO153414
SERALON	3/0	2	1 x 0,35	24	-	VO200614
SERALON	3/0	2	1 x 0,50	24	RO201714	VO201714
SERALON	3/0	2	1 x 0,75	24	RO203414	VO203414
SERALON		2,5	1 x 0,75	24	-	VO253414
SERALON	2/0	3	1 x 0,50	24	-	VO301714
SERALON	2/0	3	1 x 0,75	24	RO303414	VO303414
SERALON	0	3,5	1 x 0,50	24	-	VO351714
SERALON	0	3,5	1 x 0,75	24	-	VO353414
SERAPID	5/0	1	1 x 0,45	36	30101314	-
SERAPID	4/0	1,5	1 x 0,45	36	30151314	-
SERAPID	4/0	1,5	1 x 0,70	36	30153214	-
SERAPID	3/0	2	1 x 0,45	36	30201314	-
SERAPID	3/0	2	1 x 0,70	36	30203214	-
SERAPID	2/0	3	1 x 0,45	36	30301314	-
SERAPREN	5/0	1	1 x 0,50	24	-	CO101714
SERAPREN	4/0	1,5	1 x 0,50	24	-	CO151714
SERAPREN	4/0	1,5	1 x 0,75	24	-	CO153414

Continued on next page!

Needle designation

4 x DS-18
Multipack

DS-21

DS-25

Material	Size/USP	Size/EP	Length of suture/metre	Units per box	Code number undyed	Code number dyed
SERAPREN	4/0	1,5	1 x 0,90	24	-	CO154414
SERAPREN	3/0	2	1 x 0,50	24	-	CO201714
SERASYNTH	4/0	1,5	1 x 0,45	24	-	90151314
SERASYNTH	4/0	1,5	1 x 0,70	24	-	90153214
SERASYNTH	3/0	2	1 x 0,45	24	-	90201314
SERASYNTH	2/0	3	1 x 0,45	24	-	90301314
SUPRAMID	5/0	1	1 x 0,50	24	-	TO101714
SUPRAMID	4/0	1,5	1 x 0,50	24	-	TO151714
SUPRAMID	4/0	1,5	1 x 0,75	24	-	TO153414
SUPRAMID	3/0	2	1 x 0,35	24	-	TO200614
SUPRAMID	3/0	2	1 x 0,50	24	SO201714	TO201714
SUPRAMID	3/0	2	1 x 0,75	24	-	TO203414
SUPRAMID	2/0	3	1 x 0,50	24	SO301714	TO301714
SERALON	4/0	1,5	4 x 0,50	12	-	VO152014
SERALENE	3/0	2	1 x 0,50	24	-	LO201797
SERALON	4/0	1,5	1 x 0,50	24	-	VO151797
SERALON	4/0	1,5	1 x 0,75	24	-	VO153497
SERALON	3/0	2	1 x 0,50	24	-	VO201797
SERALON	3/0	2	1 x 0,75	24	-	VO203497
SERALON	2/0	3	1 x 0,50	24	-	VO301797
SERALON	2/0	3	1 x 0,75	24	-	VO303497
SERAPREN	3/0	2	1 x 0,75	24	-	CO203497
SERAPREN	2/0	3	1 x 0,75	24	-	CO303497
SERAFAST	4/0	1,5	1 x 0,70	24	40153215	-
SERAFAST	3/0	2	1 x 0,70	24	40203215	JO203215
SERAFAST	2/0	3	1 x 0,70	24	40303215	-
SERAFIT	4/0	1,5	1 x 0,45	36	50151315	-
SERAFIT	3/0	2	1 x 0,70	36	50203215	60203215
SERAFIT	2/0	3	1 x 0,70	36	50303215	60303215
SERAFIT	2/0	3	1 x 0,90	36	-	60304415
SERAFIT	0	3,5	1 x 0,90	36	-	60354415
SERALENE	4/0	1,5	1 x 0,50	24	-	LO151715
SERALENE	3/0	2	1 x 0,50	24	-	LO201715
SERALENE	3/0	2	1 x 0,75	24	-	LO203415
SERALENE	2/0	3	1 x 0,50	24	-	LO301715
SERALENE	2/0	3	1 x 0,75	24	-	LO303415
SERALENE	0	3,5	1 x 0,75	24	-	LO353415
SERALON	5/0	1	1 x 0,50	24	-	VO101715
SERALON	5/0	1	1 x 0,75	24	-	VO103415
SERALON	4/0	1,5	1 x 0,50	24	-	VO151715
SERALON	4/0	1,5	1 x 0,75	24	-	VO153415
SERALON	3/0	2	1 x 0,50	24	-	VO201715
SERALON	3/0	2	1 x 0,75	24	-	VO203415
SERALON	3/0	2	1 x 0,90	24	-	VO204415
SERALON		2,5	1 x 0,50	24	-	VO251715

Continued on next page!

Needle designation

Material	Size/USP	Size/EP	Length of suture/metre	Units per box	Code number undyed	Code number dyed
SERALON		2,5	1 x 0,75	24	-	VO253415
SERALON	2/0	3	1 x 0,50	24	-	VO301715
SERALON	2/0	3	1 x 0,75	24	-	VO303415
SERALON	0	3,5	1 x 0,50	24	-	VO351715
SERALON	0	3,5	1 x 0,75	24	-	VO353415
SERALON	1	4	1 x 0,50	24	-	VO401715
SERALON	1	4	1 x 0,75	24	-	VO403415
SERAPID	4/0	1,5	1 x 0,70	36	30153215	-
SERAPID	3/0	2	1 x 0,45	36	30201315	-
SERAPID	3/0	2	1 x 0,70	36	30203215	-
SERAPID	2/0	3	1 x 0,45	36	30301315	-
SERAPID	2/0	3	1 x 0,70	36	30303215	-
SERAPREN	4/0	1,5	1 x 0,50	24	-	CO151715
SERAPREN	4/0	1,5	1 x 0,75	24	-	CO153415
SERAPREN	3/0	2	1 x 0,50	24	-	CO201715
SERAPREN	3/0	2	1 x 0,75	24	-	CO203415
SERAPREN	2/0	3	1 x 0,50	24	-	CO301715
SERAPREN	2/0	3	1 x 0,75	24	-	CO303415
SERASYNTH	3/0	2	1 x 0,70	24	-	90203215
SERASYNTH	2/0	3	1 x 0,70	24	-	90303215
SERASYNTH	0	3,5	1 x 0,45	24	-	90351315
SERASYNTH	0	3,5	1 x 0,70	24	-	90353215
SUPRAMID	4/0	1,5	1 x 0,50	24	-	TO151715
SUPRAMID	3/0	2	1 x 0,50	24	-	TO201715
SUPRAMID	3/0	2	1 x 0,75	24	-	TO203415
SUPRAMID	2/0	3	1 x 0,50	24	-	TO301715
SUPRAMID	2/0	3	1 x 0,75	24	SO303415	TO303415
SUPRAMID	0	3,5	1 x 0,50	24	-	TO351715
SUPRAMID	0	3,5	1 x 0,75	24	SO353415	TO353415
SUPRAMID	1	4	1 x 0,50	24	-	TO401715
SUPRAMID	1	4	1 x 0,75	24	SO403415	TO403415
SERAFIT	2/0	3	1 x 0,70	36	-	60303216
SERAFIT	0	3,5	1 x 0,70	36	50353216	60353216
SERALENE	0	3,5	1 x 0,50	24	-	LO351716
SERALON	4/0	1,5	1 x 0,50	24	-	VO151716
SERALON	3/0	2	1 x 0,50	24	-	VO201716
SERALON	3/0	2	1 x 0,75	24	-	VO203416
SERALON		2,5	1 x 0,75	24	-	VO253416
SERALON	2/0	3	1 x 0,50	24	-	VO301716
SERALON	2/0	3	1 x 0,75	24	-	VO303416
SERALON	0	3,5	1 x 0,50	24	-	VO351716
SERALON	0	3,5	1 x 0,75	24	-	VO353416
SERALON	0	3,5	1 x 0,90	24	-	VO354416
SERAPID	2/0	3	1 x 0,70	36	30303216	-
SERAPID	0	3,5	1 x 0,70	36	30353216	-
SERAPREN	3/0	2	1 x 0,75	24	-	CO203416
SERAPREN	2/0	3	1 x 0,75	24	-	CO303416
SERAPREN	0	3,5	1 x 0,50	24	-	CO351716

DS-30



Continued on next page!

Needle designation

Material	Size/USP	Size/EP	Length of suture/metre	Units per box	Code number undyed	Code number dyed
SUPRAMID	3/0	2	1 x 0,75	24	-	TO203416
SUPRAMID	2/0	3	1 x 0,50	24	-	TO301716
SUPRAMID	2/0	3	1 x 0,75	24	-	TO303416
SUPRAMID	0	3,5	1 x 0,50	24	-	TO351716
SUPRAMID	0	3,5	1 x 0,75	24	-	TO353416
SERALON	3/0	2	4 x 0,50	12	-	VO202016
SUPRAMID	0	3,5	4 x 0,50	12	-	TO352016
SERAFIT	2/0	3	1 x 0,70	36	-	60303217
SERAFIT	0	3,5	1 x 0,70	36	-	60353217
SERAPID	2/0	3	1 x 0,70	36	30303217	-
SERAPREN	2/0	3	1 x 0,75	24	-	CO303417
SERAPREN	0	3,5	1 x 0,75	24	-	CO353417
SUPRAMID	3/0	2	1 x 0,50	24	-	TO201717
SUPRAMID	2/0	3	1 x 0,50	24	-	TO301717
SUPRAMID	2/0	3	1 x 0,75	24	-	TO303417
SUPRAMID	0	3,5	1 x 0,50	24	-	TO351717
SUPRAMID	0	3,5	1 x 0,75	24	-	TO353417
SUPRAMID	1	4	1 x 0,75	24	-	TO403417
SERALON	3/0	2	1 x 0,50	24	-	VO201717
SERALON	3/0	2	1 x 0,75	24	-	VO203417
SERALON	2/0	3	1 x 0,50	24	-	VO301717
SERALON	2/0	3	1 x 0,75	24	-	VO303417
SERALON	0	3,5	1 x 0,50	24	-	VO351717
SERALON	0	3,5	1 x 0,75	24	-	VO353417
SERALON	1	4	1 x 0,50	24	-	VO401717
SERALENE	7/0	0,5	1 x 0,50	24	-	LO05171Z
SERALON	7/0	0,5	1 x 0,50	24	-	VO05171Z
SERALON	6/0	0,7	1 x 0,50	24	-	VO07171Z
SERASYNTH	6/0	0,7	1 x 0,45	24	-	9007131Z
SERAFIT S	6/0	0,7	1 x 0,50	36	5007171S	-
SERAFIT S	5/0	1	1 x 0,45	36	5010131S	-
SERALENE	6/0	0,7	1 x 0,50	24	-	LO07171S
SERALON	7/0	0,5	1 x 0,50	24	-	VO05171S
SERALON	6/0	0,7	1 x 0,50	24	-	VO07171S
SERALON	5/0	1	1 x 0,50	24	-	VO10171S
SERAPID	6/0	0,7	1 x 0,45	36	3007131S	-
SERAPID	6/0	0,7	1 x 0,70	36	3007321S	-
SERAPREN	6/0	0,7	1 x 0,50	24	-	CO07171S
SERASYNTH	6/0	0,7	1 x 0,45	24	-	9007131S
SUPRAMID	6/0	0,7	1 x 0,75	24	-	TO07341S

4 x DS-30 Multipack



DS-35



DSS-8



DSS-10



Needle designation

DSS-13

Material	Size/USP	Size/EP	Length of suture/metre	Units per box	Code number undyed	Code number dyed
SERAFAST	5/0	1	1 x 0,45	24	4010131L	-
SERAFAST	4/0	1,5	1 x 0,45	24	4015131L	-
SERAFIT	6/0	0,7	1 x 0,45	36	5007131L	6007131L
SERAFIT	6/0	0,7	1 x 0,70	36	-	6007321L
SERAFIT	5/0	1	1 x 0,45	36	5010131L	6010131L
SERAFIT	4/0	1,5	1 x 0,45	36	5015131L	6015131L
SERALENE	6/0	0,7	1 x 0,50	24	-	LO07171L
SERALENE	5/0	1	1 x 0,50	24	-	LO10171L
SERALENE	5/0	1	1 x 0,75	24	-	LO10341L
SERALENE	4/0	1,5	1 x 0,50	24	-	LO15171L
SERALON	7/0	0,5	1 x 0,50	24	-	VO05171L
SERALON	6/0	0,7	1 x 0,50	24	-	VO07171L
SERALON	6/0	0,7	1 x 0,75	24	-	VO07341L
SERALON	5/0	1	1 x 0,50	24	RO10171L	VO10171L
SERALON	5/0	1	1 x 0,75	24	-	VO10341L
SERALON	4/0	1,5	1 x 0,50	24	-	VO15171L
SERALON	4/0	1,5	1 x 0,75	24	-	VO15341L
SERALON	3/0	2	1 x 0,50	24	-	VO20171L
SERAPID	5/0	1	1 x 0,45	36	3010131L	-
SERAPID	5/0	1	1 x 0,70	36	3010321L	-
SERAPID	4/0	1,5	1 x 0,45	36	3015131L	-
SERAPID	3/0	2	1 x 0,45	36	3020131L	-
SERAPREN	6/0	0,7	1 x 0,50	24	-	CO07171L
SERAPREN	6/0	0,7	1 x 0,75	24	-	CO07341L
SERAPREN	5/0	1	1 x 0,50	24	-	CO10171L
SERAPREN	5/0	1	1 x 0,75	24	-	CO10341L
SERAPREN	4/0	1,5	1 x 0,50	24	-	CO15171L
SERASYNTH	5/0	1	1 x 0,45	24	-	9010131L
SERASYNTH	4/0	1,5	1 x 0,45	24	-	9015131L
SUPRAMID	6/0	0,7	1 x 0,50	24	-	TO07171L
SUPRAMID	5/0	1	1 x 0,50	24	-	TO10171L



DSS-15

SERAFIT	5/0	1	1 x 0,45	36	5010131T	-
SERAFIT	4/0	1,5	1 x 0,45	36	-	6015131T
SERALENE	6/0	0,7	1 x 0,50	24	-	LO07171T
SERALENE	5/0	1	1 x 0,50	24	-	LO10171T
SERALENE	4/0	1,5	1 x 0,50	24	-	LO15171T
SERALENE	3/0	2	1 x 0,50	24	-	LO20171T
SERALON	6/0	0,7	1 x 0,50	24	-	VO07171T
SERALON	6/0	0,7	1 x 0,75	24	-	VO07341T
SERALON	5/0	1	1 x 0,35	24	-	VO10061T
SERALON	5/0	1	1 x 0,50	24	-	VO10171T
SERALON	5/0	1	1 x 0,75	24	-	VO10341T
SERALON	4/0	1,5	1 x 0,50	24	RO15171T	VO15171T
SERALON	4/0	1,5	1 x 0,75	24	-	VO15341T
SERALON	3/0	2	1 x 0,50	24	RO20171T	VO20171T
SERALON	3/0	2	1 x 0,75	24	-	VO20341T

Continued on next page!



Needle designation

DSS-18

Material	Size/USP	Size/EP	Length of suture/metre	Units per box	Code number undyed	Code number dyed
SERAPID	4/0	1,5	1 x 0,45	36	3015131T	-
SERAPREN	5/0	1	1 x 0,50	24	-	CO10171T
SERAPREN	4/0	1,5	1 x 0,50	24	-	CO15171T
SERAPREN	3/0	2	1 x 0,50	24	-	CO20171T
SERASYNTH	5/0	1	1 x 0,45	24	-	9010131T
SUPRAMID	5/0	1	1 x 0,50	24	-	TO10171T
SERAFAST	4/0	1,5	1 x 0,45	24	4015131M	-
SERAFAST	4/0	1,5	1 x 0,70	24	4015321M	-
SERAFAST	3/0	2	1 x 0,45	24	4020131M	-
SERAFIT	4/0	1,5	1 x 0,45	36	5015131M	6015131M
SERAFIT	3/0	2	1 x 0,45	36	5020131M	6020131M
SERAFIT	3/0	2	1 x 0,70	36	5020321M	-
SERALENE	5/0	1	1 x 0,50	24	-	LO10171M
SERALENE	4/0	1,5	1 x 0,50	24	-	LO15171M
SERALENE	4/0	1,5	1 x 0,75	24	-	LO15341M
SERALENE	3/0	2	1 x 0,50	24	-	LO20171M
SERALENE	3/0	2	1 x 0,75	24	-	LO20341M
SERALENE	2/0	3	1 x 0,50	24	-	LO30171M
SERALENE	2/0	3	1 x 0,75	24	-	LO30341M
SERALON	5/0	1	1 x 0,35	24	-	VO10061M
SERALON	5/0	1	1 x 0,50	24	-	VO10171M
SERALON	5/0	1	1 x 0,75	24	-	VO10341M
SERALON	4/0	1,5	1 x 0,50	24	-	VO15171M
SERALON	4/0	1,5	1 x 0,75	24	-	VO15341M
SERALON	3/0	2	1 x 0,50	24	-	VO20171M
SERALON	3/0	2	1 x 0,75	24	-	VO20341M
SERALON	2/0	3	1 x 0,50	24	-	VO30171M
SERALON	2/0	3	1 x 0,75	24	-	VO30341M
SERAPID	4/0	1,5	1 x 0,45	36	3015131M	-
SERAPID	4/0	1,5	1 x 0,70	36	3015321M	-
SERAPID	3/0	2	1 x 0,45	36	3020131M	-
SERAPID	3/0	2	1 x 0,70	36	3020321M	-
SERAPREN	4/0	1,5	1 x 0,50	24	-	CO15171M
SERAPREN	4/0	1,5	1 x 0,75	24	-	CO15341M
SERAPREN	3/0	2	1 x 0,50	24	-	CO20171M
SERAPREN	3/0	2	1 x 0,75	24	-	CO20341M
SERASYNTH	4/0	1,5	1 x 0,45	24	-	9015131M
SERASYNTH	3/0	2	1 x 0,45	24	-	9020131M
SERASYNTH	2/0	3	1 x 0,45	24	-	9030131M
SUPRAMID	4/0	1,5	1 x 0,50	24	-	TO15171M
SUPRAMID	3/0	2	1 x 0,50	24	-	TO20171M



Needle designation	Material	Size/USP	Size/EP	Length of suture/metre	Units per box	Code number undyed	Code number dyed	
DSS-21	SERALON	4/0	1,5	1 x 0,50	24	-	VO15172N	
	SERALON	4/0	1,5	1 x 0,75	24	-	VO15342N	
	SERALON	3/0	2	1 x 0,50	24	-	VO20172N	
	SERALON	3/0	2	1 x 0,75	24	-	VO20342N	
DSS-24	SERAFAST	4/0	1,5	1 x 0,45	24	4015131N	-	
	SERAFAST	3/0	2	1 x 0,45	24	4020131N	-	
	SERAFAST	3/0	2	1 x 0,70	24	4020321N	-	
	SERAFIT	3/0	2	1 x 0,45	36	5020131N	-	
	SERAFIT	2/0	3	1 x 0,45	36	5030131N	-	
	SERALENE	3/0	2	1 x 0,50	24	-	LO20171N	
	SERALENE	3/0	2	1 x 0,75	24	-	LO20341N	
	SERALENE	2/0	3	1 x 0,50	24	-	LO30171N	
	SERALENE	2/0	3	1 x 0,75	24	-	LO30341N	
	SERALON	4/0	1,5	1 x 0,50	24	-	VO15171N	
	SERALON	4/0	1,5	1 x 0,75	24	-	VO15341N	
	SERALON	3/0	2	1 x 0,50	24	-	VO20171N	
	SERALON	3/0	2	1 x 0,75	24	-	VO20341N	
	SERALON		2,5	1 x 0,50	24	-	VO25171N	
	SERALON	2/0	3	1 x 0,50	24	-	VO30171N	
	SERALON	2/0	3	1 x 0,75	24	-	VO30341N	
SERALON	0	3,5	1 x 0,50	24	-	VO35171N		
SERALON	0	3,5	1 x 0,75	24	-	VO35341N		
SERAPID	4/0	1,5	1 x 0,70	36	3015321N	-		
SERAPID	3/0	2	1 x 0,70	36	3020321N	-		
SERAPREN	3/0	2	1 x 0,50	24	-	CO20171N		
SERAPREN	3/0	2	1 x 0,75	24	-	CO20341N		
SERAPREN	2/0	3	1 x 0,50	24	-	CO30171N		
SERAPREN	2/0	3	1 x 0,75	24	-	CO30341N		
SUPRAMID	2/0	3	1 x 0,50	24	-	TO30171N		
DSS-30	SERALENE	2/0	3	1 x 0,75	24	-	LO30341P	
	SERALON	3/0	2	1 x 0,50	24	-	VO20171P	
	SERALON	3/0	2	1 x 0,75	24	-	VO20341P	
	SERALON	2/0	3	1 x 0,50	24	-	VO30171P	
	SERALON	2/0	3	1 x 0,75	24	-	VO30341P	
	SERALON	0	3,5	1 x 0,50	24	-	VO35171P	
	SERALON	0	3,5	1 x 0,75	24	-	VO35341P	
	SERAPREN	0	3,5	1 x 0,75	24	-	CO35341P	
	DSS-40	SERALON	2/0	3	1 x 0,75	24	-	VO30341Q

Further combinations on request!


Needle designation	Material	Size/USP	Size/EP	Length of suture/metre	Units per box	Code number undyed	Code number dyed
Intracutaneous suture incutifix®							
DS-15	SERALON incutifix	4/0	1,5	1 x 0,50	12	VW151713	
	SERALON incutifix	3/0	2	1 x 0,50	12	VW201713	
DS-18	SERALENE incutifix	4/0	1,5	1 x 0,50	12	LW151714	
	SERALENE incutifix	3/0	2	1 x 0,50	12	LW201714	
	SERALENE incutifix	2/0	3	1 x 0,50	12	LW301714	
DS-25	SERALON incutifix	4/0	1,5	1 x 0,50	12	VW151714	
	SERALON incutifix	3/0	2	1 x 0,50	12	VW201714	
	SERAPREN incutifix	3/0	2	1 x 0,50	12	CW201714	
DS-30	SERALON incutifix	3/0	2	1 x 0,50	12	VW201715	
	SERALON incutifix	3/0	2	1 x 0,75	12	VW203415	
	SERALON incutifix	2/0	3	1 x 0,50	12	VW301715	
	SERALON incutifix	2/0	3	1 x 0,75	12	VW303415	
SERAPREN incutifix	2/0	3	1 x 0,50	12	CW301715		
DS-30	SERALON incutifix	3/0	2	1 x 0,50	12	VW201716	
DSS-13	SERALON incutifix	5/0	1	1 x 0,50	12	VW10171L	
DSS-18	SERALENE incutifix	3/0	2	1 x 0,50	12	LW20171M	
	SERALENE incutifix	2/0	3	1 x 0,50	12	LW30171M	
DSS-18	SERALON incutifix	3/0	2	1 x 0,50	12	VW20171M	
	SERAPREN incutifix	4/0	1,5	1 x 0,50	12	CW15171M	
	SERAPREN incutifix	3/0	2	1 x 0,50	12	CW20171M	
	SERAPREN incutifix	3/0	2	1 x 0,75	12	CW20341M	
DSS-24	SERALENE incutifix	3/0	2	1 x 0,50	12	LW20171N	
	SERALON incutifix	3/0	2	1 x 0,50	12	VW20171N	
DSS-24	SERALON incutifix	2/0	3	1 x 0,50	12	VW30171N	
	SERAPREN incutifix	3/0	2	1 x 0,50	12	CW20171N	
	SERAPREN incutifix	3/0	2	1 x 0,75	12	CW20341N	
SERAPREN incutifix	2/0	3	1 x 0,75	12	CW30341N		
GS-45	SERALON incutifix	2/0	3	1 x 0,50	12	VW301735	

Further combinations on request!

PRE-CUT SUTURE THREADS
WITHOUT NEEDLE/ABSORBABLE

PRODUCT RANGE

Material	Size/USP	Size/EP	Length of suture/metre	Units per box	Code number undyed	Code number dyed
SERAFIT	5/0	1	3 x 0,45	24	5N101200	6N101200
SERAFIT	4/0	1,5	3 x 0,45	24	5N151200	6N151200
SERAFIT	3/0	2	3 x 0,45	24	5N201200	6N201200
SERAFIT	2/0	3	3 x 0,45	24	5N301200	6N301200
SERAFIT	0	3,5	3 x 0,45	24	5N351200	6N351200
SERAFIT	1	4	3 x 0,45	24	-	6N401200
SERAFIT	4/0	1,5	6 x 0,45	24	5N151500	6N151500
SERAFIT	3/0	2	6 x 0,45	24	5N201500	6N201500
SERAFIT	2/0	3	6 x 0,45	24	5N301500	6N301500
SERAFIT	0	3,5	6 x 0,45	24	5N351500	6N351500
SERAFIT	1	4	6 x 0,45	24	5N401500	6N401500
SERAFIT	3/0	2	12 x 0,45	24	-	6N207400
SERAFIT	2/0	3	12 x 0,45	24	-	6N307400
SERAFIT	0	3,5	12 x 0,45	24	-	6N357400
SERAFIT	1	4	12 x 0,45	24	-	6N407400
SERAFIT	4/0	1,5	2 x 0,70	12	-	6N157200
SERAFIT	3/0	2	2 x 0,70	12	-	6N207200
SERAFIT	2/0	3	2 x 0,70	12	-	6N307200
SERAFIT	0	3,5	2 x 0,70	12	-	6N357200
SERAFIT	1	4	2 x 0,70	12	-	6N407200
SERAFIT	4/0	1,5	5 x 0,70	12	-	6N157500
SERAFIT	3/0	2	5 x 0,70	12	-	6N207500
SERAFIT	2/0	3	5 x 0,70	12	-	6N307500
SERAFIT	0	3,5	5 x 0,70	12	-	6N357500
SERAFIT	1	4	5 x 0,70	12	-	6N407500
SERAPID	4/0	1,5	3 x 0,45	24	3N151200	-
SERAPID	3/0	2	3 x 0,45	24	3N201200	-
SERAPID	2/0	3	3 x 0,45	24	3N301200	-
SERAPID	0	3,5	3 x 0,45	24	3N351200	-
SERAPID	4/0	1,5	6 x 0,45	24	3N151500	-
SERAPID	3/0	2	6 x 0,45	24	3N201500	-
SERAPID	2/0	3	6 x 0,45	24	3N301500	-
SERAPID	0	3,5	6 x 0,45	24	3N351500	-
SERAPID	1	4	6 x 0,45	24	3N401500	-
SERAPID	2/0	3	2 x 0,70	12	3N307200	-
SERAPID	0	3,5	2 x 0,70	12	3N357200	-
SERAPID	1	4	2 x 0,70	12	3N407200	-
SERAPID	1	4	5 x 0,70	12	3N407500	-



Material	Size/USP	Size/EP	Length of suture/metre	Units per box	Code number undyed	Code number dyed
SERALON	5/0	1	3 x 0,50	24	-	VN101900
SERALON	4/0	1,5	3 x 0,50	24	-	VN151900
SERALON	3/0	2	3 x 0,50	24	-	VN201900
SERALON	2/0	3	3 x 0,50	24	-	VN301900
SERALON	4/0	1,5	5 x 0,50	24	-	VN152100
SERALON	3/0	2	5 x 0,50	24	-	VN202100
SERALON	2/0	3	5 x 0,50	12	-	VN302100
SERALON	4/0	1,5	10 x 0,50	24	-	VN152300
SERALON	3/0	2	10 x 0,50	24	-	VN202300
SERALON	1	4	2 x 0,75	24	-	VN403500
SUPRAMID	3/0	2	3 x 0,50	24	-	TN201900
SUPRAMID	2/0	3	3 x 0,50	24	SN301900	TN301900
SUPRAMID	0	3,5	3 x 0,50	24	-	TN351900
SUPRAMID	1	4	3 x 0,50	12	SN401900	-
SUPRAMID	4/0	1,5	5 x 0,50	24	-	TN152100
SUPRAMID	3/0	2	5 x 0,50	24	-	TN202100
SUPRAMID	2/0	3	5 x 0,50	12	SN302100	-
SUPRAMID	0	3,5	5 x 0,50	12	SN352100	-
SUPRAMID	1	4	5 x 0,50	12	-	TN402100
SUPRAMID	3/0	2	10 x 0,50	24	-	TN202300



Further combinations on request!

SERALON	5/0	1	1 x 100	1	-	VD108000
SERALON	4/0	1,5	1 x 100	1	-	VD158000
SERALON	3/0	2	1 x 100	1	-	VD208000
SERALON		2,5	1 x 100	1	-	VD258000
SERALON	2/0	3	1 x 100	1	-	VD308000
SERALON	0	3,5	1 x 100	1	-	VD358000
SERALON	1	4	1 x 100	1	-	VD408000
SUPRAMID	5/0 ¹	1	1 x 100	1	-	TC108000
SUPRAMID	4/0	1,5	1 x 100	1	-	TC158000
SUPRAMID	3/0	2	1 x 100	1	-	TC208000
SUPRAMID	2/0	3	1 x 100	1	-	TC308000
SUPRAMID	0	3,5	1 x 100	1	-	TC358000
SUPRAMID	1	4	1 x 50	1	-	TC408400



Further combinations on request!


PRE-CUT SUTURE THREADS
WITHOUT NEEDLE/NON-ABSORBABLE


PRODUCT RANGE

SUTURE THREADS WITHOUT NEEDLE
NON-ABSORBABLE/CASSETTES





SERAG-WIESSNER KG
Zum Kugelfang 8 - 12
95119 Naila/Germany


 + 49 9282 937-0


 + 49 9282 937-9369

Export Department:

 + 49 9282 937-230

 + 49 9282 937-9785

 info@serag-wiessner.de

 www.serag-wiessner.de